



SuperPole™

An incredibly solid, expansion fit, floor to ceiling pole. Helpful for sitting and standing. Install in virtually any room.

Features and Benefits

• Maintain Confidence and Independence.

Place support exactly where it's needed for valuable assistance beside a bed, bath, toilet or chair.

• Quick, Easy, Secure Installation.

Simply turn the screw at the bottom of the pole and the pole will expand between the floor and ceiling, creating unparalleled sturdiness and reliability.

• Versatile.

Need to remove or relocate the pole? Not a problem! Rubber pads on the floor and ceiling plates prevent marking, even on stipple ceilings.

• Innovative.

Should mobility become a greater challenge, the SuperPole™ supports a complete range of options: SuperBar™, SuperTrapeze™, and Safety Island™.



Jackscrew Installation



www.orcahealthcare.com
1-888-761-9451

Technical Specifications:



SuperPole™	90" / 229 cm x 1.5" / 3.81 cm diameter. White powder coated steel.
Grip Length	20" / 51 cm. Washable, latex free grip.
Floor to Ceiling Range	93-99" / 236-253 cm (standard).
Expansion Mechanism	1-1/8" fine cut thread. Zinc plated steel jackscrew.
Top Plate	20" / 51 cm length. Top plate rubber will not mark ceiling.
Base Plate	5" / 13 cm diameter. Non-marking, base plate rubber ribbed for better gripping.
Weight Capacity	300 lbs / 136 kg. Bariatric weight capacity 450 lbs / 205 kg.

SuperPole™ Accessories (each sold separately):

Uni-Fit Extender™	Increase floor to ceiling range from 100-120" / 254-305 cm. Not to be used with the SuperTrapeze™, Bariatric, & Ultra Models.
SuperPole™ Ultra	Longer pole accommodates floor to ceiling range from 121-140" / 307-356 cm. Not to be used with the SuperBar™, Safety Island™ or SuperTrapeze™.
Ceiling Plate Extender™	Allows standard ceiling plate 20" / 51 cm length to span 30" / 76 cm with holes for ceiling attachment to 24" / 61 cm ceiling beam centers.