Why choose any other FIXED bathlift?



AquaJoy Bathlifts Ltd Consett Business Park Consett **County Durham** DH8 6BN **United Kingdom**

email: info@aqua-joy.com www.aqua-joy.com













AquaJoy Saver & Junior Saver Bathlifts











So why choose any other FIXED bathlift?

SUITABLE FOR ADULTS AND CHILDREN

Adaptable accessories for users with varying ages and needs

BARIATRIC WEIGHT LIMIT AS STANDARD 160KG (24 STONE)

· Suitable for heavier clients

FITS ALMOST ANY STYLE OF BATH

- Ideal for deep baths maximum seat height 465mm
- Ideal for shallow baths minimum seat height 78mm

MAXIMISES LEGROOM IN THE BATH

- · Bathlift seat locates as far back as possible in bath
- · Short footprint and seat

INTERCHANGEABLE SPARE PARTS AND ACCESSORIES

Saving time, effort, stockholding and money!

Moving and Handling - made easy!

HEAVIEST PART ONLY 4.65KGS

- Reduces the risk of moving and handling injuries
- · Reduces back strain when fitting and removing

COMPACT DESIGN WHICH CURVES AROUND THE BODY

- Follows guidelines for the safe lifting of loads
- Allows an unobstructed view

EASY RELEASE SUCTION FEET WITH LIFTING BAR AS STANDARD

- Limits stooping and stretching when lifting the Bathlift in and out of the bath
- No struggling under the bathlift seat to detach individual suction feet













Easy to clean and maintain

CAN BE USED WITH OR WITHOUT COVERMATS

• Detachable and can be washed in a standard washing machine up to 60 degrees protecting against unwanted germs and bacteria

SMOOTH, CREVICE FREE SURFACE

- Dirt and used water do not get trapped
- · Makes cleaning quick and easy

STREAMLINED FRAMEWORK

• Provides open access to the Bathlift and bath for cleaning when in situ

Simple and safe to operate

LIGHTWEIGHT, FULLY WATERPROOF, FLOATING HANDCONTROL

- Comfortable to hold and operate
- Suction pads on reverse to stick to the side of the bath or to tiles for ease of use

CONVEX AND CONCAVE BUTTONS

- Soft to touch
- · Easy to operate for users with limited dexterity or visual impairment

SIMPLE CHARGING PROCESS

Pronged connectors for ease of use

BATTERY INDICATOR LIGHT

· Indicates when recharging is required

CHARGER INDICATOR LIGHT

• Indicates charging is in progress

SAFETY FEATURE

• The Bathlift will not lower unless there is sufficient charge to raise the client back to the top of the bath







